



Product: Bovine Anti-HSA (Human Serum Albumin), Affinity Purified.

Catalog Number: BHSA

Lot Number

Size: 0.5mg in 0.5ml PBS containing 0.02% sodium azide as preservative.

Immunogen: Purified human serum albumin.

Specificity: The antibody reacts with human serum albumin, and may have some cross-reactivity with albumin of other species

Suggested Dilution Range: 1:2,000-1:10,000 for ELISA, Western Blot, color substrates
1:10,000-1:40,000 for Western Blot, chemiluminescence

Short Term Storage: 2-8°C as an undiluted liquid for several weeks. After dilution, do not use for more than one day.

Long Term Storage: -20°C for up to one year.

Note: For in vitro research use only, not for diagnostic or therapeutic use.

Western blot analysis of BHSA (Bovine Anti-HSA, 1:5000 dilution) on human plasma (0.5ug) and HSA (0.1ug) shows that human serum albumin is detected by the antibody.

