



**Product:** Bovine Anti-HSA (Human Serum Albumin) Biotin Conjugate, Affinity Purified.

**Catalog Number:** BHSAB

**Lot Number**

**Size:** 0.5mg in 0.5ml PBS containing 0.02% sodium azide as preservative.

**Immunogen:** Purified human serum albumin.

**Specificity:** The antibody reacts with human serum albumin, and may have some cross-reactivity with albumin of other species

**Suggested Dilution Range:** 1:2,000-1:10,000 for ELISA, Western Blot, color substrates  
1:10,000-1:40,000 for Western Blot, chemiluminescence

**Short Term Storage:** 2-8°C as an undiluted liquid for several weeks. After dilution, do not use for more than one day.

**Long Term Storage:** -20°C for up to one year.

**Note:** For in vitro research use only, not for diagnostic or therapeutic use.

Western blot analysis of BHSAB (Bovine Anti-HSA Biotin, 1:5000 dilution) on human plasma (0.5ug) and HSA (0.1ug) shows that human serum albumin is detected by the antibody.

