



Product: Chicken Anti-GPR177, G Protein-Coupled Receptor 177, affinity purified.
Human, mouse and rat specific.

Catalog Number: OA0023

Lot Number: 80719

Size: 0.1mg in 0.1ml PBS with 0.02% Na Azide

Immunogen: Human GPR177 C-terminus peptide.

Specification: OA0023 binds to a strong band at about 55KDa on HEK293 cell line, the apparent MW may be bigger or smaller. The product may cross react with GPR177 from other mammalian species.

Suggested Usage: The product can be used at 1:2,000 to 1:5,000 for immunoblot testing. User needs to optimize the concentrations in assay. The product can be used in combination with GeneTel's anti-IgY conjugates such as RCYHRP (anti-IgY HRP) and RCYAP (anti-IgY AP) Western Blot analysis.

Short Term Storage: 2-8°C for one week.

Long Term Storage: -20°C as an undiluted liquid for up to two year.
Please feel free to call us if you have any question regarding to this product or any other GeneTel's product.



MW 293

Western Blot analysis of OA0023 chicken anti-GPR177 (1:2000) on HEK293 cell transfected with GPR177-GFP. The 293 cell protein lysate contains both transfected GPR177-GFP fusion protein (75KDa) and endogenous GPR177 protein (~55KDa).