



**GeneTel Laboratories, LLC** *Leader in Proteomics Abs and Abs Arrays*  
1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

**Product:** Rabbit Anti-Bovine beta Casein, affinity purified Abs, react with bovine milk beta casein.

**Catalog Number:** OA0108

**Lot Number:** 71209

**Size:** 0.1mg in 0.1ml PBS

**Immunogen:** Bovine beta casein peptides

**Suggested Usage:** The product can be used at 1:5,000 to 1:20,000 for both ELISA and immunoblot testing. Do not use dry milk in blocking or dilution buffer in your assay. Please use BSA (3-5% in blocking and 1-3% in dilution buffers) in your assay. User needs to try out the titration for immunohistochemical testing.

**Short Term Storage:** 2-8°C for one week.

**Long Term Storage:** -20°C as an undiluted liquid for up to two year.

**Note:** Please use BSA at 3-5% in blocking buffer and 1-3% in dilution buffer in your assay, do not use dry milk in blocking or dilution buffer.

Please feel free to call us if you have any question regarding to this product or any other GeneTel's product.



1 2 3

Western Blot analysis of OA0108 lot 71208 (lane 1) and 71209 (lane 2) at 1:10,000 dilution against 0.5ul bovine milk. Lane 3 is MW marker.