

GeneTel Laboratories, LLC

Leader in Proteomics Abs and Abs Arrays

1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

Product: Chicken Anti-ADRB3, adrenergic β 3 receptor; affinity purified; mouse and rat specific.

Catalog Number: OB0013

Lot Number: 61117

Size: 0.1ml PBS with 0.02% Na Azide, 0.5mg/ml IgY Antibodies.

Immunogen: Mouse ADRB3 peptide (LNRFDGYEGARPFPT).


Specification: The antibody reacts specifically with ADRB3 protein (about 44KD) of mouse and rat and may cross-react with ADRB3 protein of other mammalian species.

Suggested Usage: The product can be used at 1:2000-1:5000 dilutions for Western Blot analysis and for other immuno-assays. Customers need to optimize the dilution when using the Abs in IHC or ICC.

Short Term Storage: 2-8°C for one week.

Long Term Storage: -20°C as an undiluted liquid for up to two year.

Note: For in vitro research use only, not for diagnostic or therapeutic use.



Western Blot analysis of OB0013 Chicken Anti-ADRB3 Abs (1:2000) with mouse brain extract (40ug).