



GeneTel Laboratories, LLC

Leader in Proteomics Abs and Abs Arrays

1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

Product: Chicken Anti-Mouse C-Reactive Protein, affinity purified.

Catalog Number: OB0014

Lot Number: 61218

Size: 0.1ml PBS with 0.02% Na Azide, 0.5mg/ml IgY Antibodies.

Immunogen: Mouse C-Reactive protein C-terminal peptide.

Specification: The antibody reacts specifically with mouse C-reactive protein, and may have some cross-reaction with C-reactive protein from rat, human and other mammalian species.

Suggested Usage: The product can be used at 1:2000-1:5000 for both ELISA and Western Blot analysis. Customers need to optimize the dilution when using the Abs in IHC or ICC.

Short Term Storage: 2-8°C for one week.

Long Term Storage: -20°C as an undiluted liquid for up to two year.

Note: For in vitro research use only, not for diagnostic or therapeutic use.



Western Blot analysis of OB0014 Chicken Anti-mouse C-reactive protein antibodies (1:5000) on mouse plasma (2ul).