



GeneTel Laboratories, LLC *Leader in Proteomics Abs and Abs Arrays*

1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

Product: Chicken Anti-Rat C-Reactive Protein, affinity purified.

Catalog Number: OB0015

Lot Number: 61219

Size: 0.1ml PBS with 0.02% Na Azide, 0.5mg/ml IgY Antibodies.

Immunogen: Rat C-Reactive protein peptide (aa108-121).

Specification: The antibody reacts specifically with rat C-reactive protein, and may cross-react with C-reactive protein from mouse, human and other mammalian species.

Suggested Usage: The product can be used at 1:2000-1:5000 for both ELISA and Western Blot analysis. Customers need to optimize the dilution when using the Abs in IHC or ICC.

Short Term Storage: 2-8°C for one week.

Long Term Storage: -20°C as an undiluted liquid for up to two year.

Note: For in vitro research use only, not for diagnostic or therapeutic use.



Western Blot analysis of OB0015 Chicken Anti-rat C-reactive protein antibodies (1:5000) on rat plasma (2ul).