

**GeneTel Laboratories, LLC** *Leader in Proteomics Abs and Abs Arrays*  
1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

**Product:** Chicken Anti-Adiponectin; affinity purified; human, mouse and rat specific.

**Catalog Number:** OB0021

**Lot Number:** 70131

**Size:** 0.1ml PBS with 0.02% Na Azide, 0.5mg/ml IgY Antibodies.

**Immunogen:** Mouse and human Adiponectine peptide (NDSTFTGFLLYHDTN).


**Specification:** The antibody reacts specifically with Adiponectin protein (about 30KD) of human, mouse, rat and may cross-react with Adiponectin protein of other mammalian species.

**Suggested Usage:** The product can be used at 1:2000-1:5000 dilutions for Western Blot analysis and for other immuno-assays. Customers need to optimize the dilution when using the Abs in IHC or ICC.

**Short Term Storage:** 2-8°C for one week.

**Long Term Storage:** -20°C as an undiluted liquid for up to two year.

**Note:** For in vitro research use only, not for diagnostic or therapeutic use.



Western Blot analysis of OB0021 Chicken Anti-Adiponectin Abs (1:2000) with mouse plasma (1ul).