



Product: Chicken Anti-GNB3 (G-Protein $\beta 3$); affinity purified; human, mouse and rat specific.

Catalog Number: OB0029

Lot Number: 70425

Size: 0.1ml PBS with 0.02% Na Azide, 0.5mg/ml IgY Antibodies.

Immunogen: Human GNB3 peptide (GTCRQTFTGHESDIN).

Specification: The antibody reacts specifically with GNB3 protein (about 38KD) of human, mouse, rat and may cross-react with GNB3 protein of other mammalian species.

Suggested Usage: The product can be used at 1:1000-1:5000 dilutions for Western Blot analysis and for other immuno-assays. Customers need to optimize the dilution when using the Abs in IHC or ICC.

Short Term Storage: 2-8°C for one week.

Long Term Storage: -20°C as an undiluted liquid for up to two year.

Note: For in vitro research use only, not for diagnostic or therapeutic use.



Western Blot analysis with OB0029 Chicken Anti-GNB3 Abs (1:2000) on human adipose protein extract (30ug).