



Genetel Laboratories, LLC

Leader in Proteomics Abs and Abs Arrays

1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

Product: Chicken anti-phospho-MEK1/MEK2 (pS²²²/pS²²⁶), affinity purified.

Catalog Number: PMEK1

Lot Number

Size: 50ul, in PBS with 0.02% Na Azide.

Immunogen: Peptide containing phosphorylated motif (LIDSMANpSFVGTRSY) of MEK1 or MEK2 Kinase

Suggested Usage: The product can be used at 1:1,000 to 1:4,000 for immunoblot testing. Customer needs to optimize dilutions for IHC or ICC. The product can be used in combination with GeneTel's conjugates such as DRAP, DRCY3, DRFITC or GRHRP.

Short Term Storage: 2-8°C for one week.

Long Term Storage: -20°C as an undiluted liquid for up to two year.

Note: For in vitro research use only, not for diagnostic or therapeutic use.

Please feel free to call us if you have any question regarding to this product or any other GeneTel's product.

MW

50-

36-



Detection of phosphorylated MEK1/MEK2 in NGF-treated PC12 cell extract (40ug).