



Genetel Laboratories, LLC

Leader in Proteomics Abs and Abs Arrays

1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

Product: Rabbit Anti-Chicken IgY, Fc specific, Affinity Purified.

Catalog Number: RCYFc

Lot Number: 11007

Size: 0.5mg in 0.5ml PBS buffer containing 0.02% Na Azide.

Suggested Dilution Range: 1:1,000-1:5,000 for ELISA, Western Blot, color substrates
1:5,000-1:10,000 for Western Blot chemiluminescence

Preservative: 0.02% sodium azide.

Short Term Storage: 2-8°C as an undiluted liquid for several weeks. After dilution, do not use for more than one day.

Long Term Storage: -20°C for up to one year.

Specificity: The antibody reacts with constant regions of chicken IgY.

Note: For in vitro research use only, not for diagnostic or therapeutic use.
IgY is the original designation of chicken IgG like immunoglobulin.