



Genetel Laboratories, LLC

Leader in Proteomics Abs and Abs Arrays

1202 Ann Street, Madison, WI 53713 Phone 877-248-4316 Fax 608-257-1385 Email info@genetel-lab.com

Product: Chicken Anti-Ratt IgG (H&L), Affinity Purified.

Catalog Number: YRTG

Lot Number

Size: 0.5mg in 0.5ml PBS containing 0.02% sodium azide as preservative.

Immunogen: Purified Rat IgG.

Specificity: The antibody reacts with both heavy and light chains of rat IgG, and may have some cross-reactivity with IgG of other species

Suggested Dilution Range: 1:2,000-1:10,000 for ELISA, Western Blot, color substrates
1:10,000-1:40,000 for Western Blot, chemiluminescence

Short Term Storage: 2-8°C as an undiluted liquid for several weeks. After dilution, do not use for more than one day.

Long Term Storage: -20°C for up to one year.

Note: For in vitro research use only, not for diagnostic or therapeutic use.

Western blot analysis of YRTG (Chicken IgY Anti-Rat IgG, 1:4000 dilution) on rat plasma (5ug) shows that rat IgG heavy chain and light chain are detected by the antibody.



MW Rat Plasma
Marker 5ug